

THE Graft™ 

THE Graft™  Collagen

OpenTex® 

OpenTex®-TR 

Biotex® 

BioCover™ 



812, 27, Dunchon-daero 457beon-gil,
Jungwon-gu, Seongnam-si,
Gyeonggi-do, 13219, Korea
Tel: +82 2 548 1875
Made in Korea
www.purgo-biologics.com

EU Importer
Purgo Biologics Europe

1 Square Felix Bloch,
Pôle Activ' Ocean,
85300 Challans, France
Tel: +32 (0)2 28 10 61 02
E-mail: europe@purgobiologics.com
www.purgo-europe.com

EU Authorized Representatives
OBELIS S.A

Bd. Général Wahis, 53
1030 Brussels, Belgium
Tel: +32 2 732 59 54
Fax: +32 2 732 60 03
E-mail: mail@bbefe.net
www.obelis.net

THE Graft™, OpenTex®,
OpenTex®-TR, Biotex®,
THE Graft™ Collagen.

Copyright 2021. Purgo Biologics Inc.; all rights reserved . PN: EU_CA_P_20210923






“Feel the Clinical Freedom on Science and Safety”

Founded in 1999, Purgo Biologics strives to become one of the leading global companies in oral health care with its focus on safe biomaterials for soft tissue and bone regeneration.

Based on the specialized experiences accumulated by our outstanding research personnel, Purgo Research and Development Center based in Seoul is thriving to become the best in the world, specifically in the expertise of oral biomaterials for soft tissue and bone regeneration. All members in Research and Development Center are pursuing the optimized technical developments with various clinical studies, cooperative research with the governments, clinicians and educational institutions.

The solutions manufactured by Purgo are gaining fame throughout the world and Purgo’s solutions are widely accepted by global dentists from more than 30 countries.

Our production site is complying with the most international quality standards and regularly inspected by international agencies. Each production stage of our biologics solutions are controlled from the selection of the raw material to the final product.



*We had a desire.
A desire to provide Valuable & Worthwhile products for our family. That’s why we are here to let them smile shine and brightly again.*

Purgo Biologics

THE Graft™

Natural bone substitute
Page 4 to 9

THE Graft™ Collagen

Bone substitute with collagen
Page 10 to 13

OpenTex®

Non-resorbable PTFE membrane
Page 14 to 17

OpenTex® -TR

PTFE Titanium reinforced membrane
Page 18 to 21

Botex®

Suture for dental implant surgery
Page 22 to 25

BioCover™

Resorbable collagen membrane
Page 26 to 27





Science Speaks THE Graft™

THE Graft™

THE Graft™ is a natural, porcine bone mineral matrix for bone reconstruction. It is produced by removal of all organic components from porcine bone. Due to its native structure THE Graft™ likens the physical and chemical aspects of mineralized matrix of human bone. When packed into a bone defect, THE Graft™ gradually resorbs and replaced with bone during the healing process. It is available in cancellous granules packaged in vial. THE Graft™ is sterilized using gamma irradiation.

Unique proprietary manufacturing process removes very effectively potential immunogenic organic elements keeping the natural structure of the matrix.

THE Graft™ quality and safety have been scientifically demonstrated with in-vitro, in-vivo studies, large case study reports and international randomized clinical research. Systematic review and meta-analysis are conducted on THE Graft™ worldwide. [1-2]

THE Graft™ has established its fame throughout the world, both scientifically and clinically, becoming a popular bone regeneration material.



Specifications

Description	Item N°	Size / Volume ~	
THE Graft (S*) Granule 0.25 - 1.00mm	BG-A25	0.60cc	0.25g
	TG-AS25	0.25cc	
	BG-A05	1.20cc	0.50g
	TG-AS05	0.50cc	
	BG-A10	2.40cc	1.00g
	TG-AS10	1.00cc	
THE Graft (L*) Granule 1.00 - 2.00mm	BG-A20	4.80cc	2.00g
	TG-BS25	0.25cc	
	BG-B05	1.80cc	0.50g
	TG-BS05	0.50cc	
	BG-B10	3.60cc	1.00g
	TG-BS10	1.00cc	

* S : small / L : large

Indications

BONE REPLACEMENT MATERIALS	GR/CC	Extraction socket with intact socket	Extraction socket with defective socket	Minor bone augmentation	Major bone augmentation	Sinur floor elevation	Peri-implantitis
THE Graft™ Granules 0,25-1mm	0.25g~0,6cc	•	•	•			•
THE Graft™ Granules 0,25-1mm	0.50g~1,2cc	•	•	•			•
THE Graft™ Granules 0,25-1mm	1.00g~2,4cc	•	•	•	•	•	•
THE Graft™ Granules 1-2mm	0.50g~1,8cc				•	•	
THE Graft™ Granules 1-2mm	1.00g~3,6cc				•	•	



[1] Alveolar ridge regeneration of damaged extraction sockets using deproteinized porcine versus bovine bone minerals: A randomized clinical trial. 100 patients Clin Implant Dent Relat Res 2018 Jul 27. Epub 2018 Jul 27.

[2] Randomized clinical trial of maxillary sinus grafting using deproteinized porcine and bovine bone mineral. 16 Patients clin implant dent relai Res. 2017; 19(1): 140-150

«Safety and purity are an important concern when using a biomaterial»

THE Graft™ Purity [3-4-5]

Is THE Graft™ safety material?

Proprietary virus inactivation process technology.

Thanks to highly efficient manufacturing process, THE Graft™ is free from any organic components that might be potential causes of infection or immune reaction. In addition the unique process helps preserve the physical properties of THE Graft™ with its native osseous structure. A large surface area is a key requirement for graft materials, and not only results in a larger surface region available for osteoblast cells attachment but also facilitates the exchange of nutrients and waste products, it allows greater amounts of blood, proteins, and growth factors to be absorbed onto the scaffold.

THE Graft™ has a high purity.

The analysis result minimal residual protein, soft tissue, and organic bone matrix, proves that THE Graft™ is deproteinized enough for safe use.

Other than THE Graft™, such low values for organic residues are only found with bone graft material treated at high temperatures which may cause the detriment of the natural bone structure.



These results show that organic substances, including collagen and other organic compounds, were successfully removed from THE Graft™, which is thus not affected by issues associated with organic content. [3]




Is porcine bone safer than bovine?



THE Graft™ demonstrated a protein content lower than that of the natural bovine bone graft material. Bovine cancellous bone is Not Free of Zoonoses, such as BSE-Bovine Spongiform Encephalopathy. Porcine bone has a relatively low risk of zoonosis.

Less residual organic content for High purity



High purity means low organic matters

- ▶ High Surface Energy
- ▶ High Hydrophilicity



[3] Physicochemical characterization of porcine bone-derived grafting material and comparison with bovine xenografts for dental applications. Jung Heon Lee, Gyu Sung Yi, Jin Woong Lee, Deug Jeong Kim, School of Advanced Materials Science and Engineering, Sungkyunkwan University, Suwon, Korea 2SKKU Advanced Institute of Nanotechnology, Sungkyunkwan University, Suwon, Korea

[4] Process Development of a virally-safe dental xenograft material from porcine bones, Dong-Myon Kim, Ho-Chang Kang, Hyung-Joon Cha, Jung Eun Bae, and In Seop Kim, Korean Journal of Microbiology [2016] Vol. 52, No.2, pp. 140-147



[5] Effect of the calcination temperature on the composition and microstructure of hydroxyapatite derived from human and animal bone, M. Figueiredo, A. Fernando, G. Martins, J. Freitas, F.Judas, H. Figueiredo - Ceramics International Volume 36, Issue 8, December 2010, Pages 2383-2393 [2016]140-147

THE Graft™ Biocompatibility [3-4-5]

« Getting closer to human bone »


THE Graft™ is structurally similar to human bone. It has high possible level of porosity combined with a natural interconnectivity.

Safe & Biocompatible

The combination of porcine origin with the high level of purity enables predictable bone growth without risking an immunogenic reaction. In an In-Vitro study THE Graft™ was shown to encourage cell adhesion to the same extent as the compared DBBM (Deproteinized bovine bone matrix), and therefore offering optimal conditions for vital cell growth.




Porosity is an important factor in determining tissue-implant material integration. High porosity leads to a quicker absorption of liquids and cells spreading. THE Graft™ provides the optimized bone architecture for cells adhesions and tissue regeneration.



Human bone

THE Graft™

Comparison of bone structure and composition from human and animal origin




THE Graft™ High Porosity [3-4-5-6]

High porosity and early remodelling improve clinical performance.

The high porosity of THE Graft™ means a quicker absorption of fluids (e.g; blood) in comparison with DBBM. This not only facilitates the application of the material but also leads to early remodelling and improved clinical performance.

High level of porosity was demonstrated with particle pore structure test, particle size distribution test and total porosity tests.



THE Graft™ Structure :

1 -Macropores (diameter → 100µm), are necessary to form blood vessels and induce both bone growth and reorganization around the graft material.

2 -Micropores (diameter ←10 µm), are required for the penetration of body fluids, ion transportation, the attachment of osteoblasts, and the precipitation of newly formed HA.

3 -Nanopores are characterized by dimension of less than 100 nanometers pores size between grains. Nano-porosity increases bone graft permeability to the physiological fluids and cells adhesion.



Global porosity analysis :




THE Graft™ Hydrophilicity [3]

THE Graft™ consists of a unique inter-connecting pore system that ensures an efficient fluid intake and permits the migration of cells. This pore system and high surface energy enhance the osteoconduction process.

The SSA of THE Graft™ was similar with the values measured for the «bovine bone a» and significantly larger than the «bovine bone b». Considering that both THE Graft™ and «the bovine bone a» had a similar surface morphology and pore size distribution with a substantial amount of nanoscale pores. we believe that this difference in the SSA was closely related to the nano/microscale structure of the bone graft materials.


Specific surface area



THE Graft™ has shown to have higher wettability than the compared xenografts. It suggests that THE Graft™ is relatively hydrophilic and can be easily wet by body fluids after implantation. Not only protein adsorption, but also the attachment, growth, and proliferation of various types of cells, including osteoblasts, have been reported to be significantly affected by the wettability of the material surface.

This high wettability of THE Graft™ suggests that it may have advantages in terms of protein adsorption and the resulting cell adhesion and proliferation processes after implantation. The content of the organic component of THE Graft™ was somewhat lower than compared existing xenografts.

Wetting mass of the graft materials as a function of time.



This result indicates that the wettability of THE Graft™ was significantly higher than the bovine bone.



[6] Osteogenic effect of low-temperature-heated porcine bone particles in a rat calvarial defect model., Go A 1, Kim SE, Shim KM, Lee SM, Choi SH, Son JS, Kang SS . J Biomed Mater Res A. 2014 Oct;102[10]:3609-17. doi: 10.1002/jbm.a.35022. Epub 2013 Nov 18.



THE Graft™ Collagen

Biocompatible and safe natural bone grafting material, THE Graft™ Collagen is just the science itself!

THE Graft™ Collagen block is composed of porcine derived bone mineral matrix from cancellous bone and atelocollagen from porcine tendon. THE Graft™ Collagen is a bone graft intended to fill, augment, and/or reconstruct periodontal, oral, and maxillofacial defects.

THE Graft™ Collagen bone mineral matrix is similar to physical and chemical aspects of human bone mineralized matrix. Hydrated collagen components have viscosity that facilitate for blending bone mineral matrix. With this characterization, THE Graft™ Collagen can be trimmed and/or molded to the various defect shapes and can be fixed in bone defect site.

As time passes, THE Graft™ Collagen is partially transformed by the osteoclast and osteoblast cells



TCB-01



TCB-02




TCB-03

Specifications


THE Graft™ Collagen

Description	Item N°	Size/ Volume ~		
THE Graft™ Collagen Block Type	●TCB-01	0.34cc	150mg	7x7x7
	●TCB-02	0.72cc	300mg	8x9x10
	●TCB-03	1.32cc	600mg	10x11x12



THE Graft™

THE Graft™ Collagen



Collagen

1 Easily Moldable

Made of 90% of THE Graft™ granules and 10% collagen, THE Graft™ Collagen is easier to mold than THE Graft™ granules alone. Therefore, THE Graft™ Collagen has better handling property compared to that of THE Graft™, making it possible to adapt grafting materials to various shapes of the defect site with more ease.

2 Optimal osteoconductivity

While retaining better handling properties, THE Graft™ Collagen is able to form sufficient osseous tissue for implant placement and maintain natural volume and great adhesion property which lead to minimum chair time.

3 Predictable clinical results

With great hydrophilicity, THE Graft™ Collagen stabilizes the clot and aids in revascularization of the grafting material in the defect area to increase cell migration efficiency to the mineral substrate. As a result, fast bone formation can be expected, as well as a predictable clinical result.



It is advantageous for shape and space maintenance



OpenTex®

OpenTex® Non-Resorbable PTFE Membrane is a pure medical-grade polytetrafluoroethylene (PTFE) sheet with inert biological features and predictable barrier effect. Due to the smooth surface and small pore size, OpenTex® PTFE Membrane resists the incorporation of bacteria into its structure and eases the removal of the membrane.

Non-resorbable membrane is sustainable for surgical procedure with no primary closure. OpenTex® Membrane is ideal for space-making feature providing enough space for host cells to adhere to grafting materials. OpenTex® is supplied sterile for single use only and available in various sizes. [7]



The Evolution of PTFE Membrane



Non-Resorbable PTFE membrane



Indications

GBR (Guided Bone Regeneration)

- Simultaneous use of GBR membrane and implants.
- Augmentation around implant placed in immediate extraction sites or delayed extraction sockets.

GTR (Guided Tissue Regeneration)


- Filling of bone defects after root resection, removal of cysts, and removal of retained teeth.



Specifications

OpenTex®		
Item N°	Size	
OpenTex_01	24 mm x 30 mm	
OpenTex_02	17 mm x 25 mm	

Primary Closure +++



Non-Primary Closure ++OH



[7] Alveolar ridge preservation using an open membrane approach for sockets with bone deficiency: A randomized controlled clinical trial Dong-Joo Sung DOS, MSD1 | Hyun-Chang Lim DOS, PhD2 | Dong-Woon Lee DOS, PhD1 Clin Implant Dent Relat Res. 2018; 1-8

OpenTex® Main Features



Non-Resorbable

- 100% medical grade PTFE membrane.
- Biologically inert and chemically non-reactive.
- Healing procedure is not interfered with membrane absorption.



Microporous

- Promote the gingival tissue attachment.
- Enhances ease in the interstitial fluid circulation.
- Resist the bacteria infection and fibroblast cells migration.



Minimally Invasive

- Rapid recovery of soft tissue.
- Primary Closure is not necessary.
- Virtually impervious to bacteria.
- Minimum flap reflection or dissection. Safe from bacteria infection, even in the event of the exposure.



Withstands Exposure

- Protect the tissue regeneration site.
- Regenerated underlying tissue can be evaluated.
- Provide a proper environment for the growth of blood vessel and osteogenic cells.

OpenTex® Benefits



Soft Tissue Obtaining



Aesthetic Implant Restoration



Natural Saliva Passage



Minimally Invasive

OpenTex® Strengths

1 Stability :

Non-resorbable PTFE Membrane offers enough healing time to bone regenerative process.

2 Biologically inert :

PTFE is soft tissue friendly so it is ideal material as a barrier for bone regenerative process.

3 Withstands to exposure :

PTFE Membrane withstands to exposure since it is impervious to bacteria due to their barrier function.

Characteristics of OpenTex® [8]

Impervious to Bacteria

Most of Oral Bacteria is larger than 1um. OpenTex® is micro-porous material that has the pore size small enough to prevent bacterial infiltration.

Biocompatible, OpenTex® facilitates cell adhesion on the surfaces.

Test performed shows that the surface of OpenTex® is not toxic causing cells to adhere well on the surface.




General Oral Bacteria



The matter is **PORE SIZE**

24 Hours for five cells adhesion cases on OpenTex® surface

(SEM : Scanning Electron Microscope)





OpenTex® -TR

Membrane is composed of 100% polytetrafluoroethylene (PTFE) sheet and grade 1 titanium frame, which are biologically inert and tissue compatible.

OpenTex®-TR Non Resorbable PTFE Membrane with titanium frame is designed to have a suitable surface structure and porosity to prevent integration and passage of bacteria within the interstices of the material, while maintaining space for host cells adhesion to the device.

OpenTex®-TR provides a favorable environment for neovascularization and healing of defects, through repopulating the bone derived cells and protecting the bony defects from migration of the gingival tissue derived cells.


Since the adequate space maintenance is critical to this procedure, the membrane is sufficiently stiff to prevent spontaneous collapse, but also flexible enough to easily conform to tissue contours and reduce perforations of overlying soft tissue. [9]



Indications

- 01. Extraction socket reconstruction
- 02. Bone regeneration
- 03. Where primary closure isn't possible

Primary Closure



Specifications

OpenTex® -TR

Item N°	Size
OpenTex-TR_01	17 mm x 25 mm
OpenTex-TR_02	24 mm x 30 mm
OpenTex-TR_03	17 mm x 25 mm
OpenTex-TR_05	12 mm x 24 mm
OpenTex-TR_06	14 mm x 24 mm
OpenTex-TR_07	30 mm x 40 mm



[9] Clinical evaluation of vertical ridge augmentation using titanium reinforced PTFE membrane. Department of periodontology, Dental hospital service medical center. Seoul. Korea. The Korean Academy of Oral Maxillofacial Implantology Vol.22. No.1. 2018

OpenTex®-TR Main Features



Non-Resorbable



Minimally Invasive



Optimal rigidity for space maintenance



17 x 25 mm

OpenTex-TR_01



24 x 30 mm

OpenTex-TR_02



17 x 25 mm

OpenTex-TR_03



12 x 24 mm

OpenTex-TR_05



14 x 24 mm

OpenTex-TR_06



30 x 40 mm

OpenTex-TR_07

OpenTex®-TR Benefits

1 Optimal rigidity and strength for space making
OpenTex®-TR is optimal product which is able to be trimmed easily and it is solid enough for space making since it is reinforced with titanium frame.

2 Diverse embedded titanium frame
OpenTex®-TR is designed in various shapes to meet surgeon's demand.

3 Excellent tissue interaction
Its micro porous structure helps the tissue interaction.

4 Easy of use
OpenTex®-TR can be trimmed easily and also removed easily.

Characteristics of OpenTex®-TR

PTFE sheet

Grade 1 Titanium

PTFE sheet



Barrier function bacteria resistance

Predictable hard tissue integration and bone fill

- ✓ Membrane can be molded and shaped for tenting and space maintenance.
- ✓ The rigidity of the membrane is enhanced to be used for space maintenance.
- ✓ Provides additional stability in large, non-space-making osseous defects.
- ✓ Provide with little memory of Titanium frame, which enables easy placement of the membrane.
- ✓ Ability to withstand exposure.

Minimal memory, No tangle, and Superior handling



Biotex® Non-Resorbable PTFE Suture is comprised of a single-arm, non-resorbable monofilament suture with a stainless-steel surgical needle connected to the suture. The suture is uncoated, undyed and sterile for single use only, composed of 100% PTFE.

- ✓ SOFT HANDLING
- ✓ BIOLOGICALLY INERT
- ✓ NO TANGLE
- ✓ EASY KNOTTING



Indications

- Bone grafting procedures
- Periodontal surgery
- Guided tissue regeneration
- Ridge augmentation
- Implant surgery
- Soft tissue grafts

Specifications

Biotex®					
Item N°	USP Size	Length (cm)	Needle Length (mm)	Circle	Point Type
BT3019	3-0	45	19	3/8	▽
BT3016	3-0	45	16	3/8	▽
BT4016	4-0	45	16	3/8	▽
BT4019	4-0	45	19	3/8	▽
BT5016	5-0	45	16	3/8	▽
BT4013	4-0	45	13	3/8	▽
BT5013	5-0	45	13	3/8	▽
BTP4013	4-0	45	13	1/2	⊙

Biotex® Main Features & Benefits [10]

Suture

1 High pliability (PTFE)

- Tying and bending more at ease with less unintended loosening.


2 No room for little plaque

- It dispels the possibility of any bacterial infection as well as the plaque formation and any other factors that prevent healing process.

Needle


1 Slim reverse cutting needle tip

- Precision slim cut triangular needle for small penetration area and smooth suturing.
- Minimize damage to surrounding soft tissue.




2 Strong Attachment

- Advanced technology for strong needle attachment.
- Smooth and firm connection between needle and thread.
- Rapid healing process due to the reduced bleeding from needle insertion.



3 Strong Needle

- 33% higher strengths are required to bend needle in same degree compared to other product.
- High rigidity of the needle resists to bent stress during suturing.



Needle holding clip

Designed to hold the needle in place, also allows for secure and easy release of the suture needle from its package.



Tab

Allows surgeon to easily grasp and remove the suture needle from its needle holder clip.

'Race Track' shape

Designed to prevent suture from entangling, and allows easy release of the suture.

Sturdy & Flexible Transparent Cover

Protect and give clear visibility of suture and needle. Soft and sturdy cover effectively protect the suture.



Benefits

- Soft and comfortable for patients
- Soft texture for patient comfort
- Reliable closure period
- Superior handling: provides flexibility in the positioning of a square knot. Easy to tie - Easy to remove
- Nonwicking: Elimination of bacterial wicking usually associated to monofilament
- Maintains tensile strength
- PFOA free



Adaptable Resorbable Collagen Membrane



BioCover™ is a resorbable collagen membrane consisting porcine tissues which are similar to human collagen phylogenetically. BioCover™ resorbable collagen membrane offers excellent handling, easy adaptation to bone graft materials and less time consumption in surgery.


- ✓ FLEXIBLE & ADAPTABLE
- ✓ STRONG ENOUGH FOR SUTURE
- ✓ CROSSLINKED FOR DESIRED BARRIER DURABILITY



Specifications

BioCover

Item N°	Unit Size
PG0701EZC1525	15 x 25
PG0702EZC2030	20 x 30
PG0703EZC3040	30 x 40



30 x 40 mm



20 x 30 mm



15 x 25 mm



Benefits

- Biocompatible and safe
- Excellent Handling
- Great tissue adhesion
- Cell occlusive
- Strong enough to suture

Indications

BioCover™ is intended for use in periodontal and dental surgery procedures as a material for placement in the area of periodontal defect, dental implant, bone defect or ridge reconstruction to aid in wound healing post surgery. Considering BioCover™ indications and resorption time, it is recommended to combine the membrane with bone graft to new bone healing by osteoconduction (THE Graft™).

